Best-Chair® evolutionized and developed in japan!

Best-Chair developed in Japan for multi purposed use as Wheelchair and in Emergency, as Stairway Evacuation Chair. Once, disaster takes place, it is also used as Emergency Evacuation Chair safely and quickly.

*Best-Chair Improved user-friendliness, we adopted Horizontal elliptical pipes

Consideration to those who are riding

*Arm rest. Foot rest

Designed with security for rider.



*Seat

Seat height can be adjusted by pulling the seat strap at the back of the seat.



Pursue Safety and Comfortability

*Seat

Adopted a superior cushioning and resilience, comfortability.



*Lower Limb Fixing Belt

Equipped with Arm rest, Foot rest & Fixing Belt for leg.



Consideration for the operator

*Horizontal elliptical pipes

Designed with Horizontal elliptical pipes for easy grip



*Push Bar

Add to Push Bar in handle to improve operation as Wheel chair.





*Four-wheel travel

Easy turn in landing by four-wheel





Why do people choose our Best-Chair?

Best-Chair adapte Many users'opinion!

- 1 · Handle can be adjusted according to the height of operator. Can be operated by person of any
- 2 Designed with Horizontal elliptiocal pipes for easy grip.
- 3 Use in daily Wheelchair and in Emergency, as Stairway Evaccuation Chair.

Ride

*To move from bed side to either side of the best-chair.









Carry







V belt bar

*Use the V belt to slide down stair.



Operation

*Use the Kick bar before going down the stairs.









Accessories/Option





A big pouch for keeping





Dust cover with zipper

Spec/Accessories/Option

Measurements: W495xH1295xD1020 (in operation) W495xH1030XD2020 (at store) Durable load: 180kg Weight: 15.3kg Applicable standard stairway angle: 30-40 degree Material: aluminum (pipe)/ PVC (sheet) with Fire Retardant Material

Content: Best-Chair, Dust cover, Stand, bag, User Manual DVD, Abdomen & Leg fixing Band(free size), User's Manual(Warranty)

Dedicated Stand (assembling goods), hexagon socket head cap screw 2pc, dedicated tool (Hex wrench), dedicated cabinet (Special orders)

*Please read User's Manual carefully. Please operate with multiple people for those who are not used to it. Also Consider Operator's Physical strength. Depends on the situation, it is not everyone can operate.

Specifications are subject to change without notice.